

Sika AG Report No. SI 1008-533

Cleanroom® Suitable Materials

Sika-ComfortFloor

Particle (vs. PA6): GMP A Biol. Resistance: very good

Flooring & Coating

OUPLICATE DUPLICATE

Statement of Qualification





Statement of Qualification

Customer:

Material tested:

Tests performed (in accordance with CSM procedures):

Test parameters:

Test results / Classification:

Sika Deutschland GmbH Kornwestheimerstraße 103-107 70439 Stuttgart Germany

Sikafloor ComfortFloor

- Wearing course: Sikafloor-330
- Seal coat: Sikafloor-505 W RAL 7035

Detailed test information regarding the sample (serial number, color, batch number, etc.), environment and parameters used, can be obtained from the CSM test report number SI 1008-533 issued by the Fraunhofer-Gesellschaft..

- 1) Measurement of particle emission (airborne) from material when subjected
- 2) Measurement of the biological resistance
- 1) Reel-on-disc test vs. PA6; normal force 300 N.
- 2) Fungi test (Procedure A) using a suspension of spores and bacteria test (Procedure C) using a suspension of bacteria with the test specimens according to ISO 846.
- Incubation at 24°C and 95% relative humidity and visual analysis after four (4) weeks
- 1) The level of particulate contamination emitted during application of tribological stress on the material pairing specified lies within the permissible values of the corresponding Air Cleanliness Classes in accordance with EC-GMP Annex 1, in operation.

| Material pairing | Air Cleanliness Class |
|--------------------------------|--|
| Sikafloor ComfortFloor vs. PA6 | Suitable for GMP Class A (EC-GMP Annex 1) |

2) Biological resistance

| | ISO Classification (acc. to ISO 846) | CSM Classification |
|------------------------|---|--------------------|
| Fungis (Procedure A) | 0 | excellent |
| Bacteria (Procedure C) | 1 | very good |
| Overall Result | 1 | very good |

The CSM classification according to biological resistance is based on a worst-case consideration of both procedures A and C. Therefore, the numerical classification according to ISO 846 is transfered in following CSM classification:

0 = excellent3 = weak1 = very good4 = verv weak 2 = good5 = none



Standards used for the qualification:

Test environment:

Declaration:

EC-GMP Annex 1; ISO 846; VDI 2083 Part 18 Each standard stated refers to the version valid at the time of testing.

- 1) Cleanroom fulfilling Air Cleanliness Class ISO Class 1 Specifications (in accordance with ISO 14644-1); Vertical unidirectional air flow from ceiling to floor, Air flow velocity: $0.45 \,\mathrm{m/s}$; Temperature: $71.6 \,\mathrm{^\circ F} \pm 0.9 \,\mathrm{^\circ F}$, Relative humidity: 45 % ± 5 %
- 2) Microbiological S2 laboratory

The measuring equipment used for the qualification is regularly calibrated and is based on national and international standards. In the case where no national standards exist, the measuring procedure used corresponds with technical regulations and norms valid at the time of the measurement. The documents drawn up for this procedure are available for viewing.

The validity of this statement is temporary indefinite and limited to the named product. It can be verified under www.tested-device.com. For more CSM information, visit our website at www.ipa-csm.com.

Fraunhofer Institute for Manufacturing Engineering and Automation (IPA)

Department Ultraclean Technology and Micromanufacturing

Nobelstrasse 12 70569 Stuttgart Germany